



# Modify Inventory- RF Scanner

## Description

Use this task to modify inventory using the RF Scanner, before posting the count journal.

## Trigger

Perform this task when there is a need to modify inventory during cycle count.

## Prerequisites


You have completed initial cycle counts and need to make changes.

## Helpful Hints


This must be done prior to posting the journal.

## Procedure


1. Enter your User ID and press the “Enter” key.



2. Enter or scan your Password and press the “Enter” key.




3. Enter “8” for **Modify Count** and press the “Enter” key.



```
rfc Telnet Simulator
Menu Display
RF-SMART Menu
1 - Item Inquiry
2 - SO Pick Reg
3 - TO Pick
4 - TO Ship
5 - TO Receive
6 - CaseBottle Count
7 - Item Transfer
8 - Modify Count
9 - Inv Mvmt Rcpt
10 - TO Rcv - Case
10151820 - DOLC
```

4. Enter the **Journal ID** and press the “Enter” key.



```
rfc Telnet Simulator
Menu Display
Item counting
Journal ID
```




Pressing the “F1” key will provide a list of available Journal ID’s.


5. Enter or scan the item number and press the “Enter” key.




The Item unit of measure (UOM) defaults to bottles, press “F5” to change from bottle to cases and enter the case amount.



6. Enter the modified amount in the **Quantity** field.



7. When you receive “Success” on the **Item counting** screen, you have completed this task.



## Results

You have modified your inventory.